## SUPLASYN®

### SUPLASYN® 1-SHOT

Manufactured by Fabriqu  par Fabricado por Prodotto da Hersteller der Παράγεται από την Vyroba Vyroba Fabricado por Vervaardigd door Proizvajalec Производител Wyprodukowano przez

Mylan Institutional, Coill Lu, Inverin, Co. Galway, Ireland.   Mylan Pharma Group Limited, Ireland

Expiry Date d'expiration Fecha de caducidad Data di scadenza Verwendbar bis Ημερομηνία λήξης Pozujełn  do Expiraci  Prazo de validade Vennaldatum Uporabno do Дата окончания срока годности Termin waŹności

Batch number N mero de lote N mero de lote Codice del lotto Chargen Nummer Αριθμός Παρτίδας C sio s r ze C sio s r ze N mero de lote Charge Number Številka serije Номер партии Numer serii

Sterile by filter sterilisation St rilisation par filtration Esterilizado por filtraci n Sterilizzato per filtrazione Sterilisiert durch Sterilfiltration Αποστειρωμένο δια στερίλας διήθησης Sterilizov no filtraci  Sterilizov no filtraci n Est rifi por esterilizaci n de filtraci o Sterilisierte durch sterile filtrierung Sterilizirano s filtracij  Стерилизованое посредством фильтрации Sterylizowane filtracji 

Refer to instruction leaflet Veuillez vous r f rer   la notice explicative Consulte las instrucciones de uso Leggere attentamente le istruzioni per l'uso Gebrauchsanweisung lesen Βλέπετε φυλλάδιο οδηγιών P ed ro unin si p e t e n vord pro rounin Pozrite si pribalovnn letrk s pokynmi Consulte as instru es de utiliza  o Gelees de bijulter te raadplegen Prebrnite navodilo za uporabo См. список инсрукция Листок з'ясл њ s zlotk 

For single use only! dell'articolazione o sterilizzantes con amonio quaternario. Single use only! Product de un solo uso! Mono-use! Nur zur einmaligen Anwendung ! M no yia m a xronin J n k jednorazovnm pouzitin Len na jednoronove pouzitie! Para uma unica utiliza  o! Voor eenmalig gebruik! Samo za enkratno uporabo! Tolkno dla odnorazowego Применениe! Tylko do jednoroznowego uzytku!

Peel here Ouvrir ici Abri por aqui Aprire qui Hier  ffnen Αναλϊστε εδω Odlepte zde Odlepte tu Descole aqui Hier  ffnen  ffnen lukaj Otkryvat' zdes' Odklej tutaj

Store between 4°C and 25°C A conserver entre 4°C et 25°C Almacenar entre 4°C y 25°C Conservare tra 4°C e 25°C Packungen zwischen 4°C und 25°C aufbewahren Na φυλάσσεται σε θερμοκρασία μεταξύ 4°C και 25°C Uchov vanje pri teplot ch mezi 4°C a 25°C Uchov vanje pri teplot  4°C aŹ 25°C Conserve entre 4°C e 25°C Bewaren tussen 4°C en 25°C Shranjujte pri temperaturi med 4°C in 25°C Хранит' при температур  от +4°C до +25°C Przechowywac w temp. od 4°C do 25°C

#### English

#### SUPLASYN® SUPLASYN® 1-SHOT

STERILE SODIUM HYALURONATE SOLUTION

HYALURONIC ACID (HA) is a normal component of the synovia and plays a central role in maintaining the physiological internal environment of the joint.
COMPOSITION: Viscoselactique solution of a defined molecular weight of purified hyaluronic acid, produced by fermentation.

Each syringe contains:

	SUPLASYN	SUPLASYN 1-SHOT
Hyaluronic acid sodium salt	20 mg	60 mg
Excipients q.s	2 ml	6 ml

**PROPERTIES:** Hyaluronic acid is the prototype of a wide range of saccharide biopolymers (glycosaminoglycans), important components of all extracellular tissue structures, including cartilage and synovial fluid. The active substance of SUPLASYN is a hyaluronic acid of defined molecular chain length with a high degree of purity. The introduction of SUPLASYN into the synovial space will assist in the normalization of the joint following arthrocentesis.

**INDICATIONS:** SUPLASYN and SUPLASYN 1-SHOT are indicated for the symptomatic treatment of osteoarthritis. SUPLASYN has been shown to be beneficial in osteoarthritis for the management of pain and improvement in physical function of joints. More than one joint may be treated at the same time.

**DOSAGE AND ADMINISTRATION:** Depending upon joint size, up to 6 ml may be administered intra-articularly. **SUPLASYN:** The recommended schedule for SUPLASYN supplied in a 2ml syringe is 1 injection per week for 3 weeks, but up to 6 may be given depending on patient's condition. **SUPLASYN 1-SHOT:** Intended for single administration. Use strict aseptic technique. Discard any unused portion of the syringe. To use the pre filled syringe, remove the Luer lock cap, attach a suitable cannula (recommended is 21 – 25 G depending on joint) and secure it by turning slightly. **GRADUATION ON THE SYRINGE LABEL IS TO BE USED AS A GUIDE ONLY.**

**CONTRAINDICATIONS/PRECAUTIONS:** Do not administer to patients with known hypersensitivity reactions. Respect usual precautions and contraindications for any intra-articular injection. Do not inject intra-vascularly. **SUPLASYN 1-SHOT:** Intended for single administration. Use strict aseptic technique. Discard any unused portion of the syringe. To use the pre filled syringe, remove the Luer lock cap, attach a suitable cannula (recommended is 21 – 25 G depending on joint) and secure it by turning slightly. **GRADUATION ON THE SYRINGE LABEL IS TO BE USED AS A GUIDE ONLY.** Transient short duration pain may occur following intra-articular introduction. The affected joint may show a mild local reaction like pain, feeling of heat, hyperthermia, redness, effusion, irritation, and swelling/inflammation. If these symptoms occur, rest the affected joint and apply ice locally. Symptoms subside within days for most of the patients. In some cases, mild local reactions such as pain, irritation, swelling/joint inflammation and effusion may be significantly enhanced and much more severe as an expression of hypersensitivity. In such cases, a therapeutic intervention could be necessary, e.g. aspiration of joint fluid. Local adverse reactions could be accompanied by systemic reactions such as fever, chills, or cardiovascular reactions, and in rare cases anaphylactic reactions. In extremely rare circumstances, rash/itching, urticaria, synovitis, and a drop in blood pressure have been reported following the administration of SUPLASYN. Discontinue use if adverse reactions are experienced. Avoid using SUPLASYN with sterilising or sanitising agents containing quaternary ammonium salts solutions.

**WARNING: KEEP OUT OF REACH OF CHILDREN. DO NOT USE IF BLISTER IS DAMAGED. TO BE USED BY A PHYSICIAN ONLY. FOR SINGLE USE ONLY. DISCARD UNUSED PORTION OF SYRINGE. REUSE MAY CREATE A RISK OF CONTAMINATION AND/OR CAUSE PATIENT INFECTION OR CROSS INFECTION.**

**PACKAGING:** Available in 2ml or 6ml syringes.

**STORAGE:** Store between 4°C and 25°C. DO NOT FREEZE.

Bring to room temperature before injection.

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Datum priprave besedila: januar 2017

Datum zadnje revizije: prva izdaja

Multi Packaging Solutions		171218	A
Customer	Mylan	Dimensions (mm)	520 x 330
Product	SUPLASYN 1 SHOT	Pages	2 pp
Item Code	0953L100	Drawing	520X330/F/41
Part No.		Pharmacoce/2D/Barcode	N/A
		Colours FRONT (or over)	Black
		Colours BACK (or inner)	Black
		Operator	Yarnold, Paul
		Date	31/01/2017 15:55:12
<p><small>Attention: The colours displayed in this proof are for electronic registrational purposes only. Please refer to the PANTONE® colour guides for accurate colour matching. By opening this file, you are signifying full approval of design and use. To ensure accurate viewing and printing of this PDF file, Overprint Preview must be activated for accurate colour representation. MPS Leaflet PDF files are designed to work with Adobe Acrobat 6 or higher and are PDF/X-01a:2001 compliant.</small></p>			