SUPLASYN® SUPLASYN® 1-SHOT

Manufactured by
Fabriqué par
Fabricado por
Prodotto da
Hersteller
Παράγεται από την
Vytroba
Vyrobca
Fabricado por
Vervaardigd door
Proizvajalec
Производителъ
Wyprodukowano przez

Mylan Institutional, Coill Rua, Inverin, Co. Galway, Ireland.
© Mylan Pharma Group Limited, Ireland

Expiry
Date of expiration
Fecha de caducidad
Data di scadenza
Verwendbar bis
Ημερομηνία λήξης
Použitelné do
Expíraciá
Prazo de validade
Vervaldatum
Uprabno do
Дата окончания срока годности
Termin ważności

Batch number
Número de lot
Número de lote
Codice del lotto
Chargen Number
Αριθμός Παρτίδας
Číslo šarže
Císlo šarže
Número de lote
Charge Number
Številka serije
Номер партии
Numer serii

Sterile by filter sterilisation
Stérilisation par filtration
Esterilizado por filtración
Sterilizzato per filtrazione
Sterilisiert durch Sterilfiltration
Αποστειρωμένο δια στερεάς διήθησης
Sterilizováno filtrací
Sterilizovane filtráciou
Estéril por esterilización de filtración
Steriliseate door steriele filtering
Sterilizrano s filtracij
Стерилизация посредством фильтрации
Sterylizowany filtracją

Refer to instruction leaflet
Veuillez vous référer à la notice explicative
Consulte las instrucciones de uso
Leggere attentamente le istruzioni per l'uso
Gebruiksaanwijzing lezen
Βλέπε φύλλο οδηγιών
P ed použitím si p e t ie návod pro použití
Pozrite si príbalový leták s pokynmi
Consulte as instruções de utilização
Gelieve de bijsluiter te raadplegen
Preberite navodilo za uporabo
См. листок интрукций
Zapoznaj się z ulotką

For single use only!
A usage unique!
Producto de un sólo uso!
Monouso!
Nur zur einmaligen Anwendung!
Môno για μία χρήση!
Jen k jednorázovému použitiu!
Len na jednorazovú použitie!
Para uma única utilização!
Voor eenmalig gebruik!
Samt za enkratno uporabo!
Только для однократного применения!
Тільки до одноразового узытку!

Peel here
Ouvrir ici
Abrir por aquí
Aprire qui
Hier öffnen
Ανοίξτε εδώ
Odklepte zde
Odklepte tu
Descolle aqui
Hier orenen
Odpriete tukaj
Открывать здесь
Odklej tutaj

Store between 4°C and 25°C
A conserver entre 4°C et 25°C
Almacén entre 4°C y 25°C
Conservare tra 4°C e 25°C
Packungen zwischen 4°C und 25°C aufbewahren
No φυλάσσεται σε θερμοκρασία μεταξύ 4°C και 25°C
Uchovávaťe pri teplote medzi 4°C a 25°C
Uchovávaťe pri teplote 4°C až 25°C
Conserve entre 4°C e 25°C
Bewaren tussen 4°C en 25°C
Shranjujte pri temperaturi med 4°C in 25°C
Хранити при температурі от +4°C до +25°C
Przechowywać w temp. od 4°C do 25°C

English

SUPLASYN® SUPLASYN® 1-SHOT

STERILE SODIUM HYALURONATE SOLUTION

HYALURONIC ACID (HA) is a normal component of the synovia and plays a central role in maintaining the physiological internal environment of the joint.
COMPOSITION: Viscoclastic solution of a defined molecular weight of purified hyaluronic acid, produced by fermentation.

Each syringe contains:

	SUPLASYN	SUPLASYN 1-SHOT
Hyaluronic acid sodium salt	20 mg	60 mg
Excipients q.s	2 ml	6 ml

PROPERTIES: Hyaluronic acid is the prototype of a wide range of saccharide biopolymers (glycosaminoglycans), important components of all extracellular tissue structures, including cartilage and synovial fluid. The active substance of SUPLASYN is a hyaluronic acid of defined molecular chain length with a high degree of purity. The introduction of SUPLASYN into the synovial space will assist in the normalization of the joint following arthrocentesis.

INDICATIONS: SUPLASYN and SUPLASYN 1-SHOT are indicated for the symptomatic treatment of osteoarthritis. SUPLASYN has been shown to be beneficial in osteoarthritis for the management of pain and improvement in physical function of joints. More than one joint may be treated at the same time.

DOSAGE AND ADMINISTRATION: Depending upon joint size, up to 6 ml may be administered intra-articularly.
SUPLASYN: The recommended schedule for SUPLASYN supplied in a 2ml syringe is 1 injection per week for 3 weeks, but up to 6 may be given depending on patient's condition.
SUPLASYN 1-SHOT: Intended for single administration. Use strict aseptic technique. Discard any unused portion of the syringe. To use the pre filled syringe, remove the Luer lock cap, attach a suitable cannula (recommended is 21 – 25 G depending on joint) and secure it by turning slightly.
GRADUATION ON THE SYRINGE LABEL IS TO BE USED AS A GUIDE ONLY.

CONTRAINDICATIONS/PRECAUTIONS: Do not administer to patients with known hypersensitivity reactions. Respect usual precautions and contraindications for any intra-articular injection. Do not inject intra-vascularly. SUPLASYN should not be used in patients presenting an inflammation/irritation of the joint, since adverse events more commonly occur in patients with already existing joint inflammation/irritation. As no clinical evidence is available on the use of Hyaluronic Acid in children, pregnant and lactating women, treatment with SUPLASYN is not recommended in these patients. The patient should rest 24-48 hours after the injection and avoid any strenuous activity over the full course of the treatment.

Transient short duration pain may occur following intra-articular introduction. The affected joint may show a mild local reaction like pain, feeling of heat, hyperthermia, redness, effusion, irritation, and swelling/inflammation. If these symptoms occur, rest the affected joint and apply ice locally. Symptoms subside within days for most of the patients. In some cases, mild local reactions such as pain, irritation, swelling/joint inflammation and effusion may be significantly enhanced and much more severe as an expression of hypersensitivity. In such cases, a therapeutic intervention could be necessary, e.g. aspiration of joint fluid. Local adverse reactions could be accompanied by systemic reactions such as fever, chills, or cardiovascular reactions, and in rare cases anaphylactic reactions. In extremely rare circumstances, rash/itching, urticaria, synovitis, and a drop in blood pressure have been reported following the administration of SUPLASYN. Discontinue use if adverse reactions are experienced. Avoid using SUPLASYN with sterilising or sanitising agents containing quaternary ammonium salts solutions.

WARNING: KEEP OUT OF REACH OF CHILDREN. DO NOT USE IF BLISTER IS DAMAGED. TO BE USED BY A PHYSICIAN ONLY. FOR SINGLE USE ONLY. DISCARD UNUSED PORTION OF SYRINGE. REUSE MAY CREATE RISK OF CONTAMINATION AND/OR CAUSE PATIENT INFECTION OR CROSS INFECTION.

PACKAGING: Available in 2ml or 6ml syringes.
STORAGE: Store between 4°C and 25°C. DO NOT FREEZE.

Bring to room temperature before injection.

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