EPIPEN Jr Adrenaline (Epinephrine) 150 μg/0.3 mL Auto-Injector

EpiPen® Jr. DIRECTIONS FOR USE

Keep EpiPen® Jr. away from the face and eyes.

Carefully tilt and slide the EpiPen® Jr.

Auto-Injector out of its protective carry

Check the contents of the EpiPen® Jr.

Auto-Injector through the "viewing window". Make sure the solution is clear, colourless

and sediment-free. If the "viewing window" is obscured and the needle cover is extended,

the EpiPen® Jr. has already "fired" and cannot

Make a fist around the EpiPen® Jr.

Auto-Injector, with your thumb nearest to

the blue safety release (see FIGURE 1).

Pull off the blue safety release with your

Hold the leg firmly and place the EpiPen®

Jr. Auto-injector at a 90° angle (right angle) to the outer part of your child's

thigh. Hold the orange end about 5 cm

Push the orange end FIRMLY into the child's outer mid-thigh until a "click" is heard or felt. Administer through clothing if necessary. Hold the EpiPen® Jr. Auto-Injector firmly against the thigh for

approximately 3 seconds (see FIGURE 3).

Carefully take the EpiPen® Jr. Auto-Injector

away from the thigh. The orange needle

cover will extend to cover the needle (see

Check the "viewing window" - it should now be obscured. If it is not obscured, the EpiPen® Jr. has not "fired" or worked. If the

"viewing window" is not obscured, repeat

Note the time EpiPen® Jr. [adrenaline

Auto-Injector even though some adrenaline (epinephrine) is left inside it.

Do not replace the blue safety release. Your used EpiPen® Jr. Auto-Injector with

extended orange needle cover will not fit

Seek further medical attention for yourself immediately.

10) You CANNOT RE-USE the EpiPen® Jr.

11) Do not tamper with the orange needle cover.

12) See 'IMMEDIATELY AFTER USE' box.

(epinephrine)] was given.

away from the thigh (see FIGURE 2).

other hand (also shown in FIGURE 1)

HOW TO USE EpiPen® Jr. Auto-Injector

be used again.

Grip only the mid-section of the EpiPen® Jr. Auto-Injector. Never put thumb.

stops the EpiPen® Jr. from "firing" or working.

Always point the orange end (which holds the needle) away from any part of the body, except the planned injection site. This is usually the outer mid-thigh.

Do not "click" on either end of the EpiPen® Jr. Auto-Injector. The EpiPen® Jr. Auto-Injector does not work like a retractable ball-point pen.

Do not inject into hands, feet, ears, nose, buttocks, genital area or into a blood

vessel. If this does happen accidentally you must immediately seek medical

fingers or hand over the orange end. The orange end contains the need Do NOT remove the blue safety release until ready to use. The blue safety release

Guide to using EpiPen® Jr.

Read this page BEFORE an emergency occurs. If you have any questions about this leaflet or about EpiPen® Jr. Auto-Injector, ask your doctor or pharmacist. Keep this leaflet with the EpiPen® Jr. Auto-Injector.

IMPORTANT INFORMATION

EpiPen® Jr. is for emergency treatment of children who have a severe allergic reaction (anaphylaxis) caused by hypersensitivity to (for example):

- · stinging insects
- latex
- other allergens

· insect bites.

A severe allergic reaction, otherwise known as anaphylaxis, is life-threatening. Typical symptoms include: flushing of the face

- collapse
- officulty breathing and/or wheezing

   irregular or faint pulse
- swelling (eg face, lips, tongue or throat) 
   vomiting, diarrhoea, stomach cramps.
   skin rashes, hives, itching

Selection of the appropriate dose is determined according to body weight.

The EpiPen® Jr. Auto-Injector is intended for children who weigh from 15 kg up to 30 kg, have a severe allergic reaction, and have been assessed and advised by their doctor or pharmacist as needing the EpiPen® Jr. Auto-Injector. There may be instances where the doctor will recommend the use of EpiPen® Jr. in some adults, based on careful medical assessment.

EpiPen® Jr. is life-saving in the emergency management of severe allergic reactions. It should not be withheld from anyone who is considered to need it. Discuss with your doctor or pharmacist any precautions or potential interactions. EpiPen® Jr. Auto-Injector:

- automatically injects one x 0.3 mL pre-measured dose of adrenaline (epinephrine). This dose provides 150 µg of adrenaline (epinephrine).
- · may be self-administered or administered by a carer to the child
- · cannot be reused.

For further information see the Consumer Medicine Information (CMI) for EpiPen® Jr. Auto-Injector. The CMI can be obtained from your pharmacist or downloaded from www.mylan.com.au. For New Zealand consumers, please download the CMI from www.Medsafe.govt.nz/consumers/cmi.

Do NOT remove the blue safety release until you need, and are ready to use,

The safety release stops EpiPen® Jr. from working.

8

EpiPen® Jr. is only for the immediate EMERGENCY TREATMENT of a SEVERE ALLERGIC REACTION. Your doctor or pharmacist will explain when your child must use EpiPen® Jr.

If your child has the signs and symptoms described by your doctor or pharmacist:

- Use EpiPen® Jr. immediately, through clothing if necessary. Carefully follow 'DIRECTIONS FOR USE'.

Keep calm, and still.

If your child has been stung by an insect:

Use EpiPen® Jr. immediately, through clothing if necessary.

- If the sting is easily accessible, remove with your fingernails. Do not squeeze, pinch or push the sting deeper into the skin. Ice packs or sodium bicarbonate
- soaks may then be applied to the stung area. Keep warm, calm, and as still as possible

Care and Storage of EpiPen® Jr. Auto-Injector Keep available ready for use at all times.

- Handle carefully, so as not to make it "fire" accidentally.
- Never place thumb, fingers or hands over the orange end as this contains the
- Do not remove blue safety release until ready to use. The blue safety release stops EpiPen® Jr. from "firing" or working.

Storage

- Keep in a cool dry place below 25°C. Temperature excursions between 15°C to
- Do not place in the fridge. Do NOT leave EpiPen® Jr. in your car. Extreme temperatures may cause the auto-injector unit to malfunction
- Protect from light.
- Keep in the protective carry tube until needed.

Things you must do

- Check the expiry date on your EpiPen® Jr. Auto-Injector. Replace it BEFORE this
- Each month, check the contents through the "viewing window" of your EpiPen® Jr. Auto-Injector. Make sure the solution is clear, colourless and sediment-free. If it's not, the EpiPen® Jr. should not be used. Return it to your doctor or
- pharmacist for safe disposal. Obtain a new EpiPen® Jr. immediately.

EpiPen® Jr. may be prescribed by your doctor or purchased directly from a pharmacy. Ask your doctor or pharmacist how you can get an EpiPen® Training Device to

It is important to regularly ask your doctor or pharmacist to demonstrate use with an EpiPen® Training Device

EpiPen® Jr. is distributed in Australia by:

Alphapharm Pty Limited (ABN 93 002 359 739) Level 1, 30 The Bond, 30-34 Hickson Road Millers Point NSW 2000 Australia Phone: (02) 9298 3999

Medical Information Phone: 1800 028 365 www.mylan.com.au

and distributed in New Zealand by Mylan New Zealand Limited

Auckland New Zealand Phone: 0800 168 169

Made in USA EpiPen® is a registered trademark

of Mylan, Inc.

· difficulty breathing, shortness of breath fast or noticeable heart beat

EpiPen® Jr. may cause side effects such as:

IMMEDIATELY AFTER USE

· Contact your doctor or

Either

- skin colour: pale or red; flushing · sleeplessness, restlessness
- nausea or vomiting
- dizziness

LFJR150/6

- · anxiety, fear

for inspection and safe disposal in a rigid sharps disposal unit

- throbbing sweating
- weakness

· hallucinations

 headache coldness

if you notice these, or other worrying symptoms, talk to your doctor or pharmacist.

EpiPenJr CMI.A

3020

Call an ambulance (telephone 000 in Australia or 111 in New Zealand). Explain that your child has had EpiPen® Jr. [intramuscular adrenaline

Give your used EpiPen® Jr. Auto-Injector, to the doctor's office or to a hospital

197 00 mm

inside margins 5 mm

Finished DIMENSIONS 30mm x 50mm



## **CCL Label**

**MELBOURNE** 

120 Merrindale Drive, Croydon South, VIC 3136 Australia Phone +61 3 9751 7100 www.ccllabel.com

151871 Job No: LFJR150/6 Code:

Version: 1 Operator: AH 15/06/18 - Prepress Proof Preparation Proof Generation from supplied files

## **TERMS AND CONDITIONS**

Where the purchaser approves any artwork or proofs submitted by the company, the company will not be liable for any errors or inaccuracies subsequently discovered in the goods or any work performed or produced by the company in the course of producing the goods."

I have checked this proof for proceed with manufacturing. No changes are required.

Customer Signature: Date:

> Keyline/dieline will not print Inkjet prints are indicative only, and should be matched to relevant pantone colours for printing.

> > Artwork scale 100% - actual size

FIGURE 1  showing a person preparing an price of preparing an price of preparing and 4 show the steps as the price of prepared of the price of prepared of the price of prepared of the price of the pri	

50mm







FIGURE 3



FIGURE 4